
Process G blocks

Introduction

Here is how you process G-blocks < 1 kb from IDT.

Materials

>

>

Procedure

Processing G blocks from IDT

1. Put the g-blocks in the small centrifuge for 1 min at 3000 rpm.
2. Pre-heat the thermo-mixer to 50C. And put a thermometer in.
3. Add 20 ul of sterile water to each g block.
4. Briefly vortex and spin on the table top centrifuge (a few seconds at each).
5. Put in the thermo-mixer at 50C for 20 min. Make sure the thermometer is at 50C
6. Vortex and then tabletop centrifuge for a few sec each.