
Bleeding Mice for transduction check

Introduction

How do do a cheek bleed and check that transduced cells are in the blood.

Materials

› Materials

- › Bleeding tubes located at Vida's desk with the fluid inside of them.

Procedure

Preparing Blood

1. Get blood tubes - 1 for each mouse.
2. Label blood tubes so that you know which mouse is bled into each.
3. Bleed animal into blood tube with 1 ml of PBS + EDTA.
4.  **PAUSE** Store blood at 4 C for one dahy if you are not going to use it that day.
5. Spin the blood at 6000 rpm for 4-5 min.
6. Aspirate and discard supernatant.
7. Add 1 ml Ack and resuspend the pellet.
8. Incubate at room temp for 5 - 10 min.
9. Spin at 6000 rpm for 4-5 min.
10. Aspirate and discard supernatant. Be careful not to aspirate the pellet.
11. Resuspend pellet in 200 ul FACS buffer.
12. Transfer to a labelled 96 well V-bottom plate.
13. Make antibody master mix by preparing a 1:200 dilution of each antibody in FACS buffer.
50 ul of FACS buffer per sample + 100 ul in addition (just in case)
14. Spin 96-well plate at 2000 rpm for 2 min.

-
15. Discard supernatant by flicking once over the sink.
 16. Add 50 ul antibody master mix to the wells and resuspend. Keep on ice.
 17. Incubate on ice or at 4 C in the dark for 15 min.
 18. Add 150 ul FACS buffer.
 19. Spin at 2000 rpm for 2 min and discard the supernatant by flicking once over the sink.
 20. Add 200 ul FACS buffer and resuspend.
 21. Spin at 2000 rpm for 2 min and discard supernatant by flicking once over the sink.
 22. Add 200 ul of FACS buffer, resuspend cells, and transfer to FACS tubes on ice.