
Making LB Buffer

Introduction

How to make LB buffer

Materials

› Materials

- › LB granules
- › LB flask
- › autoclave tape
- › plastic bucket

Procedure

Procedure

1. Take 25g of the LB granules and put them in one of the big LB bottles.
2. Put autoclave tape on the side
3. Grab one of the flat plastic buckets.
4. Go to the autoclave room and check if there is an open autoclave. Sign name and check in time.
5. Fill the LB bottle with 1L water. Specifically the clean water in the autoclave room.
6. Shake the LB bottle until the granules are dissolved.
7. Fill the plastic bucket with a solid amount of water, maybe like 1/4-1/3 full.
8.  **CRITICAL UNSCREW THE CAP A BIT ON THE LB BOTTLE. Don't screw it off all the way.**
9. Place the LB bottle into the bucket with the autoclave.
10. Start program 9 on the autoclave.
11. Come back in 1.5 hours.

01:30:00



12. Leave LB partially unscrewed until it is fully cooled down and then screw on fully,