

Plating PlatE from frozen

Project: Spatial_Perturb

Author: Alexander Monell

Status: In progress

Entry Created On: 17 Apr 2025 22:20:04 UTC

Entry Last Modified: 18 May 2026 22:17:54 UTC

Export Generated On: 13 Jul 2026 18:16:41 UTC

THURSDAY, 4/17/2025

Plating Plat-E Cells from Frozen

1. Check if we have D10 media.
 2. Get a 50 ml conical, put 30 ml D10 and put in water bath.
 3. Use a 15 mL conical tube and add 10 mL D10.
 4. Put the conical tube in a 37°C water bath.
 5. Move liquid nitrogen container from the southwest rack – 4th box from top.
 6. Grab Plat-E cells.
 7. Place the frozen Plat-E vial in the water bath until thawed (1–2 min max).
 8. Add 1 mL warmed media dropwise.
 9. Transfer back into 15 mL tube.
 10. Spin – can count right here.
 11. Plate in 2 Plat-E plates,
 - 2.5×10^5 cells per 10 cm plate.
- Don't change the culture medium during the first three days. It is normal to see some cells floating after the first 24 hours.
 - Don't culture cells to complete confluency. Split cells 4X to 6X every two to three days when the culture reaches 70-90% confluency.