

MC-38 Cell Line**Jeffrey Schlom, PhD, National Cancer Institute/NIH – 2/12/15 (Revised 4/21/16)****Growth Conditions:**

- Cells are adherent
- Maintain cells at 80% confluency
- Cells are normally split on Tuesdays and Fridays
- Split cells 1:10 for use in 2-3 days
- Split 1:40 for routine Tues/Fri maintenance
- (trypsin 1X for 2-3 min is sufficient)

Medium:

Dulbecco's modified MEM with 10% fetal bovine serum, 2mM glutamine, 0.1 mM nonessential amino acids, 1 mM sodium pyruvate, 10 mM Hepes, 50ug/ml gentamycin sulfate, pen/strep.

Formulation:

500 mL DMEM (4.5g/l glucose)
Remove 78.5 ml, discard

Add:

50 mL FBS (for 10%) **Use BGS**
5 mL non-essential amino acids (10mM or 100x stock)
5 mL sodium pyruvate (100 mM stock)
5 mL Glutamine (200 mM stock)
0.5 mL Gentamicin (50 mg/mL stock)
5 mL Hepes (1 M stock)
5 mL Pen/Strep (100x stock)