

From Hedrick Lab.

VSV stock production

Reagents:

The original protocol suggests to use DM-5, we used DM-10 (DMEM with 10% FCS, + pen-strep)

Preparation

- Two days earlier prepare BHK-21 cells
 - Take a confluent mid sized flask (T75)
 - Trypsinize well
 - Seed a 1/33 and 1/66 split into a new T75 flask
 - Optionally count the cells and plate out 5×10^5 cells and 1×10^6 cells

Day of the infection:

- Pick the TC flask that is 90% confluent
- Thaw the virus stock on ice
- Consider a T75 flask to contain 5×10^6 total cells
- Aim at an MOI of 0.05 (range 0.1-0.01)
- Dilute the virus so that at the end 1ml of DM10 contains the desired amount of virus particles
- Remove the old culture media
- Add the 1ml of media containing the virus, make sure that the 1ml flows over all the cells a few times,
- Incubate at 37°C for 1 hour, tilt all 15-20 min
- Then add 5 ml of warm DM-10
- Harvest the supernatant at 16-18 hours post infection
- Spin the virus supernatant (to remove cells), aliquot 50ul and store at -80° C