

How to Use Guava for cell count

- ✓ Sign up the machine in advance.
- ✓ For the first user of the day, please **CLEAN** before use! (Click **Clean and Shut Down** from **Main Menu**, then click **Start Cleaning**, it'll take ~17min)
- ✓ Please **filter your samples** before use!
- ✓ **Always load the tubes (W1-W6) with appropriate buffer** (please follow the instructions shown on **Tray Message** dialog box!)
- ✓ For the first time non-Lu lab users, please contact **Ling-Li**(lllin@ucsd.edu, 858-822-2734) or Li-Fan before you use the machine.

Procedures:

1. Turn on Guava machine, start **CytoSoft** software
Run **Clean and Shut Down** before use!
2. Dilute samples with Guava Viacount Buffer in a 96 U-bottom plate (minimum total final volume per well is 100ul; use 200ul if your cells are in low concentration)
3. Open **WorkEdit** Software to set up your plate,
Select the wells, enter the number of **Events**(the default is 1000)
☒check **Acquire this sample**
☒check **Wash both capillary and mixer"**
☒check **Mix the samples for 3 secs at High rpm"**
Enter your **Dilution Factor** and **Total Volume** of your original sample, save file
4. Click **Guava ViaCount** from **Main Menu**
Click **Start worklist**
Load 96 well plate, check tubes (fill the tubes with appropriate solution)
Select your worklist file and save file (for data file)
Click **Adjust Settings**,
Select one of the samples to adjust the threshold, then
Click **Next Step**
Click **Resume** to start reading
5. Open **Guava Excel** software for results
6. If some is going to run the machine within the next half hour, run a quick clean. If not, run **Clean and Shut down** before turning off the machine