



Exercise 4 FoxP3 Reference Control Staining Protocol

You are planning to perform an assay for detection of FOXP3. Please indicate the steps you will follow to stain **beads** to be used as controls for surface antigens AND for intracellular antigens

Cells	Beads for Surface Antibodies	Beads for FoxP3 IC Antibody
0. Wash Beads	X	X
1. Add Surface Abs	X	
2. Incubate/Wash	X	X
3. Fix/Perm	X	X
4. Wash (with Perm Buffer)	X	X
5. Add IC Ab		X
6. Incubate/Wash	X	X
7. Acquisition	X	X